3601 Evaluation of a Novel Chromogenic Medium for the Isolation and Differentiation of Salmonella and Shigella spp.

D. Berger, R. Clasen, A. Cuna, J. Knox, J. Hardy; Hardy Diagnostics, Santa Maria, CA.

Revised Abstract

Several selective and differential media are available for the isolation of Salmonella and Shigella spp., but these media do not distinguish enteric pathogens from non-pathogenic, non-lactose-fermenting organisms. Isolation on these media require excessive colony picking to screen for possible pathogens, most of which turn out to be falsepositive normal flora, such as Proteus, Providencia, Morganella, and Citrobacter spp. While chromogenic media formulations are available for Salmonella spp., none exist for Shigella spp.

To alleviate these labor-intense workups, Hardy Diagnostics' novel chromogenic media, HardyCHROM[™] SS, allows for the selective isolation and differentiation of Salmonella and Shigella spp. from non-pathogenic enteric bacteria (both lactose and non-lactose-fermenting organisms).

The Microbial Diseases Laboratory, California Department of Public Health, inoculated a total of 80 enteric isolates using HardyCHROM[™] SS in order to evaluate the ability of this medium to accurately recover and differentiate Salmonella and Shigella spp.

HardyCHROM[™] SS yielded 100% recovery and 100% appropriate colony colors for all of the Salmonella and Shigella spp. tested. Some low H₂S-producing Salmonella strains appeared turquoise as expected. Based on these findings, HardyCHROM[™] SS, utilizing a novel, patented chromogenic substrate (under exclusive license), can be employed as a cost effective, reliable method for the selective isolation and differentiation of Salmonella and Shigella spp. from non-pathogenic enteric bacteria in stool cultures.

Introduction

Salmonella and Shigella spp. continue to be a major cause of disease and illness due to food-borne and waterborne infection.⁽¹⁻⁶⁾ Many formulations of culture media (such as HE, SS, and XLD) have been developed to isolate and differentiate Salmonella and Shigella spp. from non-pathogenic enteric bacteria.^(4,6,7) Most formulations rely on carbohydrate fermentation, pH indicators, and an indicator system for the detection of hydrogen sulfide.⁽⁷⁾ These media are made selective by the addition of bile salts and can also differentiate between Salmonella, Shigella, and lactose-fermenting organisms. However, colonies of non-lactose-fermenting organisms that are non-pathogenic can appear similar in appearance to Salmonella and Shigella and must be subjected to further testing by using Triple Sugar Iron (TSI) Agar, Lysine Iron Agar (LIA), or Kligler Iron Agar (KIA).^(4,6-9) Screening of primary plates or secondary plates inoculated from enrichment broths often requires the inoculation of large numbers of secondary screening tubes and/or the use of costly automated identification systems.

The use of chromogenic substrates (chromogens) in media formulations has increased greatly in the last several years. Chromogens, when broken down by specific bacterial enzymes, will result in colored colonies. Previously, chromogenic formulations were available for Salmonella spp., but not for Shigella spp.^(5-7,10,11)

HardyCHROM[™] SS allows for the selective isolation and differentiation of both Salmonella and Shigella spp. from non-pathogenic enteric bacteria (both lactose and non-lactose-fermenting organisms). Differentiation of Salmonella and Shigella spp. from non-pathogenic bacteria is accomplished by three mechanisms: chromogenic reactions, carbohydrate fermentation, and hydrogen sulfide production. HardyCHROM[™] SS provides better differentiation of colonies obtained from clinical samples and enrichment procedures, resulting in less secondary screening of isolates and less false-positive results.

Materials & Methods

All of the strains tested were retrieved from the Microbial Diseases Laboratory, California Department of Public Health.

- 80 total strains were tested
- 2 genera were tested (Salmonella and Shigella)
- 16 different species, serotypes, or groups were tested

Suspension Preparation and Inoculation

- Isolates were grown in tryptone or brain heart infusion broth for approximately four hours in a water bath
- A loopful of each suspension was used to streak HardyCHROM[™] SS plates for isolation
- Plates were incubated at 35 °C for 18-20 hours

Kesulfs						
Genus	Species, Serotype, or Group	Number of Isolates Tested	Number of Isolates Positive for Growth on HardyCHROM [™] SS	Number of Isolates Yielding Expected Color Reaction on HardyCHROM [™] SS	% Isolates Positive (Growth)	% Isolates Positive (Color)
Salmonella	<i>enterica</i> subsp. salamae (Group II)	3	3	3	100%	100%
Salmonella	<i>enterica</i> subsp. arizonae (Group IIIa)	3	3	3	100%	100%
Salmonella	<i>enterica</i> subsp. diarizonae (Group IIIb)	2	2	2	100%	100%
Salmonella	<i>enterica</i> subsp. houtenae (Group IV)	3	3	3	100%	100%
Salmonella	serotype Typhi	5	5	5	100%	100%
Salmonella	serotype Paratyphi A	2	2	2	100%	100%
Salmonella	serotype Heidelberg	1	1	1	100%	100%
Salmonella	spp. (serotype not determined)	19	19	19	100%	100%
Salmonella	group B	1	1	1	100%	100%
Salmonella	group C2	1	1	1	100%	100%
Salmonella	group D	1	1	1	100%	100%
Shigella	sonnei	13	13	13	100%	100%
Shigella	flexneri	12	12	12	100%	100%
Shigella	boydii	8	8	8	100%	100%
Shigella	dysenteriae	4	4	4	100%	100%
Shigella	provisional serotype	2	2	2	100%	100%
Total		80	80	80	100%	100%

Salmonella enterica (ATCC[®] 14028) colonies growing on HardyCHROM[™] SS. Incubated aerobically for 24 hours at 35 °C.

Conclusion

- turquoise or colorless for Shigella dysenteriae)

References

- in 2006. J. Clin. Microbiol.; 49:232-42.

- 40:2999-3003.
- Washington, D.C.

- J. Appl. Microbiol.; 103:2046-55.

Acknowledgement

We wish to thank the Microbial Diseases Laboratory for their participation in this study.

Shigella sonnei (ATCC[®] 9290) colonies growing on HardyCHROM[™] SS. Incubated aerobically for 24 hours at 35 °C.

• HardyCHROM[™] SS yielded 100% recovery of Salmonella and Shigella spp. tested

• 100% of Salmonella spp. yielded the expected color reactions (turquoise, with or without black centers) • 100% of Shigella spp. yielded the expected color reactions (turquoise for Shigella sonnei, flexneri, and boydii;

HARDY DIAGNOSTICS

A Culture of Service[™]

 Based on these findings, HardyCHROM[™] SS can be employed as a cost effective, reliable method for the selective isolation and differentiation of Salmonella and Shigella spp. from non-pathogenic enteric bacteria

. Preliminary FoodNet Data on the Incidence of Infection with Pathogens Transmitted Commonly Through Food—10 States, 2009. MMWR Morb Mortal Wkly Rep, 2010.; 59:418-22.

2. Xia, S., et al. 2011. Prevalence and characterization of human Shigella infections in Henan Province, China,

3. Berger C.N., et al. 2011. Salmonella enterica strains belonging to O serogroup 1,3,19 induce chlorosis and wilting of Arabidopsis thaliana leaves. *Environ. Microbiol.*; [Epub ahead of print].

4. Procop, G.W., et al. 2008. A single-tube screen for Salmonella and Shigella. Am. J. Clin. Pathol.; 130:284-9. 5. Church, D.L., et al. 2010. Clinical and economic evaluation of BBL™ CHROMagar™ Salmonella (CHROMSal) versus subculture after selenite broth enrichment to CHROMSal and Hektoen enteric agars to detect enteric Salmonella in a large regional microbiology laboratory. Diagn. Microbiol. Infect. Dis.; 68:13-9.

6. Maddocks, S., et al. 2002. Comparison of CHROMagar Salmonella medium and xylose-lysine-desoxycholate and Salmonella-Shigella agars for isolation of Salmonella strains from stool samples. J. Clin. Microbiol.;

7. van Dijk, S., et al. 2009. Evaluation and implementation of a chromogenic agar medium for Salmonella detection in stool in routine laboratory diagnostics. J. Clin. Microbiol.; 47:456-8. 8. Murray, P.R., et al. 2007. Manual of Clinical Microbiology, 9th ed. American Society for Microbiology,

9. Reller, L.B. and S. Mirrett. 1975. Motility-indole-lysine medium for presumptive identification of enteric pathogens of Enterobacteriaceae. J. Clin. Microbiol.; 2:247-52.

10. Schönenbrücher, V., et al. 2008. A comparison of standard cultural methods for the detection of foodborne Salmonella species including three new chromogenic plating media. Int. J. Food. Microbiol.; 123:61-6. 11. Perry, J.D. and A.M. Freydière, et al. 2007. The application of chromogenic media in clinical microbiology.