



Buy One Box of Sterisart Canisters Get One Box of TSB or FTM with Septum Free







Membrane Model	Description				
Sartorius Sterisart Canisters					
16469GBD	Universal Luer adapter and long dual-needle spike, sterile vented				
16469GSD	Universal Luer adapter and long dual-needle spike, sterile vented, with septum				
16475GBD	Two dual-needle spikes of different length, one is sterile vented				
16475GSD	Two dual-needle spikes of different length, one is sterile vented with septum.				
16476GBD	Short dual-needle spike, sterile vented				
16477GBD	Long needle, side opening, with solid pointed tip, non-coring, included sterile venting needle				
16478GBD	Female Luer Lock connection				

FTM 10ml	FTM 20ml	FTM 100ml
FTM U66	FTM U47	FTM U427
TSB 10ml	TSB 20ml	TSB 100ml
TSB K182	TSB U38	TSB U46

Need help finding the right free media?

Website: www.HardyDiagnostics.com/MembraneFiltration.html

Phone: 1-800-266-2222 x 5599

Hardy Membrane Filtration Mini Catalog Click here to view catalog or scan QR Code



Customer Redemption Form

To receive your free item, please fill in the following information (print please) and fax this page along with your proof of purchase to Hardy.

Name			
Title			
Company			
Address			
City	State/Prvn	Zip/PC	
		·	
Model #		Serial #	



Headquarters

1430 West McCoy Lane Santa Maria, CA 93455 800.266.2222 :phone 805.346.2760 :fax Sales@HardyDiagnostics.com www.HardyDiagnostics.com

Distribution Centers

Santa Maria, California Olympia, Washington Salt Lake City, Utah Phoenix, Arizona Dallas, Texas Springboro, Ohio Lake City, Florida Raleigh, N. Carolina Albany, New York